



Immersion thermostat with digital control "Digiterm S-150"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5°C TO 150°C.

High precision temperature control

SAFETY:

STANDARD DIN 12876. 1. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET 2. LOW LIQUID LEVEL FLOAT

FEATURES

- Maximum working volume to maintain maximum temperature: 20 litres.
- A clamp and finger screw attach the unit to the tank, alternatively the extendible support bridge accessory can be used.
- Minimum fixing depth to the tank: 14 cm.
- Rear connection for additional OT100 probe. Automatic selection.
- Cooling coil.

CONTROL PANEL

1. Set temperature display indicator.
2. Bath temperature/time display indicator.
3. "SET" button.
4. Set value decrease/increase button.
5. Start/Stop button.



MODELS	Part No.	Control range °C	Stability °C	Height / Width / Depth (exterior) cm	Pump Pressure / Flowrate H ₂ O	Power W	Weight Kg
DIGITERM S-150	3000543	Amb +5 to 150	±0,1	28 18 19	150 mbar-12 l/min.	2000	3,9

Accessories (See page. 104).