



Immersion thermostat with TFT touch screen "Digiterm TFT-200"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5°C TO 200°C.

With external refrigeration from -20 to 200.

High precision temperature control

SAFETY:

- STANDARD DIN 12876. 1. SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET
2. OVER TEMPERATURE ALARM - 3. LOW LIQUID LEVEL FLOAT

FEATURES

Maximum working volume to maintain maximum temperature: 20 litres.

A clamp and finger screw attach the unit to the tank, alternatively the extendible support bridge accessory can be used.

Minimum fixing depth to the tank: 14 cm.

RS-232, USB and Ethernet output.

Cooling coil.

CONTROL PANEL

1. Set/ current temperature indicator.
2. Pump flow indicator (%).
3. Stop button.
4. Start button.
5. Push up parameter.
6. Push down parameter.
7. Set / current time indicator in hours / minutes.
8. Return button.



MODELS	Part No.	Control range °C	Stability °C	Height / Width / Depth (exterior) cm	Pump Pressure / Flowrate H ₂ O	Power W	Weight Kg
DIGITERM TFT-200	3000544	amb.+5 to 200	200 °C ±0,1	28 18 19	150 mbar-12 l/min.	2060	4,5

Accessories (See page. 104).