



Metallic thermostat dry block “Clinic-Bloc” and “Bio-Bath 2”

CLINIC-BLOC: FIXED TEMPERATURE OF 37 °C WITH CHECK THERMOMETER.

BIO-BATH 2: TEMPERATURE RANGE: FROM AMBIENT +5 TO 100 ° C.

SAFETY: NORM EN 61010. OVER TEMPERATURE CUT-OUT FITTED.

“CLINIC-BLOC”

Capacity 20 tubes (5x4) of 75x13 mm. VAC. Part No. **7001569**

“CLINIC-BLOC” 1537

Capacity 15 cuvettes (5x3) of 10x10 mm. Part No. **7001570**

MODEL

Part No.	Height / Width /Depth (exterior) cm	Stability °C	Power W	Weight Kg
7001569	11 18 28	±0.5	10	2.4
7001570	11 18 28	±0.5	10	2.4

“BIO-BATH 2”

Capacity 8 tubes (5x4) of 75x13 mm VAC.

Adjustable temperature and time.

MODEL

Part No.	Height / Width /Depth (exterior) cm	Stability °C	Power W	Weight Kg
7001562	11 11,3 14,5	±0,3	5	0,5



For VAC tubes (75x13)
Part No. 7001569



For cuvettes (10x10)
Part No. 7001570



**NEW
DESIGN**