




# Dry Block heater for tubes "Multiplaces"

CAPACITY: THREE BLOCKS.



## ACCESSORIES

**Metal Blocks:** 55 high x 95 wide x 75 mm deep.

Part No.

- 7000306** For 30 x 6 mm  $\emptyset$  tubes (total 3 blocks = 90 tubes).
- 7000716** For 20 x 1.5ml, Eppendorf tubes. (total 3 blocks = 60 tubes.)
- 7031220** For 20 x 12 mm  $\emptyset$  tubes (total 3 blocks = 60 tubes).
- 7031612** For 12 x 16 mm  $\emptyset$  tubes (total 3 blocks = 36 tubes).
- 7003208** For 8 x 20 mm  $\emptyset$  tubes (total 3 blocks = 24 tubes).
- 7001256** For 6 x 25 mm  $\emptyset$  tubes (total 3 blocks = 18 tubes).
- 7000717** Blank block without holes, can be customised to customer requirements.
- 7001474** Block for 15 cuvettes 10 x 10. Capacity: 3 blocks.
- 7001475** Block for 96 well microplates of 1.2 ml. Capacity: 1 block.

## MODEL

Part No.	Temperature °C	Stability °C	Homogeneity %	Set error %	Resolution °C	Height / Width /Depth (exterior) cm	Power W	Weight Kg
<b>7471200</b>	30 - 200	±0.75	±1.5	±2	1	11 31 31	700	7.4

Supplied complete with thermometer 0-200 °C and extraction tool for blocks.