



Water and oil baths “Termobat”

ADJUSTABLE TEMPERATURE FROM AMBIENT+5 °C TO 110 °C. STABILITY ± 2 °C TO 200 °C. STABILITY ± 3 °C.

SAFETY:

CONFORMS TO THE DIN STANDARD 12877 BY INCORPORATING A MANUALLY RESET OVER TEMPERATURE SAFETY CUT-OUT.

APPLICATIONS

To be used with a rotary evaporator for flasks up to 10 litres or as a stand alone bains marie.

FEATURES

Tank interior made of AISI 304 stainless steel, with an epoxy coated base.

The heater element is embedded in a protective INCOLOY stainless steel shield.

CONTROL PANEL

Main power switch, with luminous indicator lamp.

A hydraulic temperature thermostat regulator is synchronised with an indicator lamp.

MODELS

Part No.	Maximum Temperature °C	Capacity litres	Ø (usable) cm	Height cm	Ø (exterior) cm	Height cm	Power W	Weight Kg
3000587	110	11	29	17	29	25	2100	5.5
3000588	200	11	29	17	29	25	2100	5.5

Large capacity



ACCESSORIES MADE OF AISI 304 STAINLESS STEEL

Reduction ring cover lid of 210, 180, 145, 115 and 80 mm Ø. With a location hole for a thermometer.

Part No. **7000589**



Heater element cover shelf.

Part No. **6000055**



Thermometer.

Part No.

1095100 up to 100 °C.

1095250 up to 250 °C.