



Bath circulators “Digiterm-S-150”

CUVETTES SET OF 20 AND 27 L. CAPACITY, WITH THERMOSTATIC EQUIPMENT INCLUDED.
FURNISHING TEMPERATURES FROM AMB+5 TO 150 °C.
STABILITY: $\pm 0,1$ °C. RESOLUTION: 0.1 °C.

NEW
DESIGN

REGULATION
SYSTEM
A

SAFETY:

SAFETY STANDARD DIN 12876.

SAFETY ADJUSTABLE OVER TEMPERATURE CUT OFF WITH MANUAL RESET.
LOW LIQUID LEVEL CUT OFF.

CONTROL PANEL

Display temperature indicator of the bath.

Set temperature indicator display.

Heater function light indicator.

Start/Stop button.

“SET” button (temperature/time).

Push button to decrease set value.

Push button to increase set value.

PT 100 external probe connection.



MODELS	Part No.	Temperature Range °C	Cuvettes Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Pressure mbar	Pump Pump rate l/min	Power W	Weight Kg
DIGITERM-S-150-20	3000548	150	20	15 30 29	36 55 35	150	12	2000	9.2
DIGITERM-S-150-27	3000549	150	27	20 30 29	41 55 35	150	12	2000	10.8

ACCESSORIES see page 104.