



## Water bath "Aquaterm"

ADJUSTABLE TEMPERATURE FROM 30 °C TO 110 °C. STABILITY  $\pm 2$  °C.

### *Economic option*

#### FEATURES

AISI 304 stainless steel tank with epoxy coated base.  
Heater element embedded in high temperature, corrosive resistant INCOLOY.

#### CONTROL PANEL

Temperature regulator thermostat switch, synchronised with an indicator lamp.

#### MODEL

Part No.	Capacity litres	Ø (usable)	Height cm	Ø (exterior)	Height cm	Power W	Weight Kg
<b>3000916</b>	1.5	20	11.5	22	20.5	640	1.3



#### ACCESSORIES IN AISI 304. STAINLESS STEEL

Perforated shelf.

Part No. **6000998**

Flat Lid with 5 concentric ring set.

Part No. **3000917**