



Boiling bath "Baheer"

ADJUSTABLE TEMPERATURE FROM 30 °C TO 110 °C. STABILITY ± 3 °C.

FEATURES

Rapid boiling (6-7 minutes).

AISI 304 stainless steel tank and epoxy coated base.

Heater element embedded in INCOLOY stainless steel.

CONTROL PANEL

Temperature regulator thermostat switch, synchronised with an indicator lamp.

MODEL

Part No.	Capacity litres	Ø (usable)	Height cm	Ø (exterior)	Height cm	Power W	Weight Kg
7000540	1	14	11	15	19	640	1



ACCESSORIES MADE OF STAINLESS STEEL

Lid with bakelite handle.

Part No. **7000572**

Rack 24 haemolysis tubes up to 13 mm diameter.

Part No. **7004002**

Rack for 24 tubes up to 18 mm diameter.

Part No. **7000541**

Perforated shelf.

Part No. **6000033**