

Manufactured O.E.M



Poupinel dry heat sterilizer "Drytime II"

ADJUSTABLE TEMPERATURES FROM 50 °C UP TO 200 °C.
STABILITY: ± 6 °C.

NEW
DESIGN

APPLICATIONS

For quick surgical sterilization of diverse instruments
surgical odontological, etc.

FEATURES

Heating by shielded elements in the base which provide
a rapid temperature rise.

Flap door.

Inner chamber in AISI 304 stainless steel.

Removable tank with extraction clamps.

Epoxy-coated outer casing.

SAFETY

Over temperature cut out incorporated. EN.61010
Standard.

CONTROL PANEL

Mains switch.

Mains indicator lamp.

Hydraulic thermostat for temperature control.

Timer 0 to 120 min. with automatic off.



MODEL

Part No.	Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
2000912	1,5	6,5 25 12	15 31 20	300	4