



Hand operated centrifuges

FEATURES

For 15 ml and butyrometer tubes.

Comes with a table holding clamp.

Adaptable heads for 2 or 4 x 15 ml tubes or for 2 butyrometer tubes.

MODELS

Part No.	Comes with head for	Max. speed r.p.m.	Max. R.c.f.	Weight Kg
7000014	2 tubes	1500	360	1.1
7000015	4 tubes	1500	360	1.2
7000016	2 butyrometers	1500	370	1.2

Comes with head and metal buckets.



ADDITIONAL ACCESSORIES

Head for two tubes.

Part No. **7000017**

Head for four tubes.

Part No. **7001388**

Bucket for 15 ml tubes.

Part No. **7001500**

Head for two butyrometers.

Part No. **7000018**