



Mixer shaker for tubes "Vortex-Vib"

SHAKE BY MIXER HEAD.

Strong and economic.

FEATURES

With stirring selector by pressure, or fixed in continuous. Starting switch, 4 interchangeable platforms for multiple tubes and 50 ml flasks. Base with 4 suction cups for full adherence.

MODEL

| Part No. | Frequency r.p.m. | Height /Width /Depth (exterior) cm | Power W | Weight Kg |
|----------------|---------------------|---------------------------------------|------------|--------------|
| 7001725 | 2800 | 18 12,5 18 | 50 | 4 |

