



Chest Freezers “Combatep”

CONTROLLABLE TEMPERATURE RANGE FROM -14 °C TO -28 °C.
 CONTROLLABLE TEMPERATURE RANGE FROM -10 °C TO -40 °C.

APPLICATIONS

Suitable for preservation and freezing of sample.

FEATURES

Epoxy covered external case. Interior of aluminium, with guide supports for baskets or trays.

The interior of the door is made of easy to clean, strong, chemically resistant, moulded ABS.

Drop down lockable door with compensating support hinges. Hermetically sealed compressor, silent and vibration free.

Automatic interior light when the door is opened. (all models except Combatep S-30.)

Digital electronic temperature control.
 N.T.C. Sensor.

CONTROL PANEL

Temperature control.

Digital live and preset temperature display.

Push-button temperature selector.



MODELS	Part No.	Capacity litres	Freezing capacity in 24hrs Kg	Temperature up to °C	Height / Width / Depth (interior) cm			Height / Width / Depth (exterior) cm			Power W	Weight Kg
COMBATEMP S-30	2101040	206	28	-28	74	63	53	92	75	71	160	38
COMBATEMP L-30	2101041	343	38	-28	74	100	53	92	113	71	200	56
COMBATEMP S-40	2101181	215	30	-40	63	89	41	92	113	76	230	57
COMBATEMP L-40	2101182	352	40	-40	65	117	50	92	137	81	400	72

Note: Part No. 2101040 is supplied with 1 basket. Part No. 2101041 is supplied with 2 baskets. Part No. 2101181 is supplied with 3 baskets and Part No. 2101182 is supplied with 2 baskets.

SPARE PARTS Baskets

For Combatep. Part No.	2101040 / 2101041	2101181	2101182
Baskets	1001810	1001811	1001812