



Upright freezers “Templow” S, M,L and YC



APPLICATIONS

Preservation and storage of products that require to be preserved at low temperatures.

Upright freezer “Templow” S

DIGITAL TEMPERATURE CONTROL,
FIXED TEMPERATURE: -15 °C, -18 °C, -21 °C, -25 °C AND -32 °C.
8 DRAWERS WITH INDIVIDUAL FREEZER COMPARTMENTS.

FEATURES

Epoxy covered metallic external case.
 Strong, molded, chemically resistant ABS interior.
 Manual defrost with water evaporation system.
 Door can be hung to open left or right.

CONTROL PANEL

Mains power ON/ OFF switch.
 Temperature thermostat regulator.
 Temperature indicator lamp.
 Audible and visual alarms.
 Push button alarm reset/stop.
 Push button with quick freeze indicator.



Upright freezer “Templow” M

TEMPERATURE RANGE FROM -14 °C TO -25 °C.
7 INDIVIDUAL FREEZER COMPARTMENTS WITH 14 BASKETS.

FEATURES

Epoxy covered external case.
 Strong, chemically resistant, ABS molded case.
 Specially coated “Rilsan”, wire mesh baskets.
 Manual defrost, with water evaporation system.

The door can be hung to open left or right.
 Lockable door.

CONTROL PANEL

Temperature thermostat regulator.
 Digital temperature display.



Upright freezer “Templow” L

WITH FAN ASSISTED AIR CIRCULATION.
TEMPERATURE RANGE FROM -10 °C UP TO -35 °C.
COMES WITH 4 SHELVES.

FEATURES

External and internal case made completely of AISI 304 stainless steel, with water sump and drain.
 Chemically resistant.
 Acclimatized hermetically sealed CFC free compressor.
 Automatic defrost with electrical water evaporation system.
 Regulable shelves.

Adjustable front level supports.
 The door can be hung to open, left or right.
 Door lock.

CONTROL PANEL

Mains power ON/ OFF switch.
 Electronic temperature regulator.
 Digital temperature display.



Refrigerated cabinet with double compressor “Templow” YC

UPPER CASE: TEMPERATURES FROM 2°C UP TO 8°C. CAPACITY: 189 LITRES.
LOWER CASE: TEMPERATURES FROM -10°C UP TO -40°C. CAPACITY: 100 LITRES.
3 SHELVES AND 1 DRAWER IN UPPER CASE AND 3 DRAWERS IN LOWER CASE.

CARACTERÍSTICAS

External case covered in epoxy.
 Stainless steel internal easy for cleaning, with also stainless steel wire basket.
 Double compressor that offers maximum stability and constancy of temperature.
 Gradual adjustment of the trays for a versatile use of the interior.
 Door with lock and key.
 USB output for data recording.

CONTROL PANEL.

Digital electronic control.
 Alarms of low and high temperature, door open and power failure.
 Backup battery for alarms.
 Upper case temperature indicator display.
 Lower case temperature indicator display.
 Parameters value increase/decrease.



MODELS	Part No.	Capacity litres	Range °C	Height / Width / Dept (interior) cm	Height / Width / Dept (exterior) cm	Power W	Weight Kg	
TEMPLOW S	2101278	257	-15 to -32	151 44 42	184 60 63	110	81	Comes with 8 drawers
TEMPLOW M	2101279	513	-14 to -25	151 61 56	186 75 75	345	109	*Comes with 14 baskets
TEMPLOW L	2101280	600	-10 to -35	155 51 65	215 70 83	600	130	*Comes with 4 shelves
TEMPLOW YC	2101288	289	**2 to -40	- - -	183 70 64	360	136	Comes with 3 shelves and 4 drawers

*Comes with steel PVC laminated shelves and baskets. **in different compartments