



Storage refrigerators “Stocklow” LS, GS and LSC

WITH FAN ASSISTED CIRCULATION.

MODEL LS, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C.

MODEL GS, WITH DOUBLE GLAZED ANTI MIST DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C.

MODEL LSC, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM -17 °C TO -22 °C.



APPLICATIONS

Preservation and storage of products that require below ambient temperatures.

FEATURES

- Chemically resistant.
- Tropical hermetic compressor.
- Automatic defrost with electrical water evaporation system.
- Support rack and shelves.
- Adjustable front level supports (Models LS and GS).
- Supplied with wheels (Model LSC).
- The door can be hung to open left or right.
- Door lock.

There are two versions:

Stocklow LS and LSC which is completely made of metal, both interior and exterior, AISI 304 stainless steel door.

Stocklow GS inside made of AISI 304 stainless steel and epoxy white external case with glazed door that enables internal view.

CONTROL PANEL

- Mains power, ON/OFF switch.
- Electronic temperature regulator.
- Digital temperature display.
- Interior light only (Double glazed door model).

SPARE PARTS

Código **1001858** Set 2 guides + additional tray



Model Stocklow LS and LSC, with metallic door.



Model Stocklow GS, with metallic door.

MODELS	Part No.	Capacity litres	Door	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Number of shelves	Power W	Weight Kg
STOCKLOW LS	2101283	600	Metallic	154 61 63	213 73 84	3	250	120
STOCKLOW GS	2101284	600	Glass	154 61 63	213 73 84	3	250	125
STOCKLOW LSC	2101289	600	Metallic	154 61 63	213 73 84	3	500	125

Comes complete with steel shelves, PVC laminated.