



Refrigerated cabinet "Medilow-XS"

FOR TEMPERATURES FROM 2°C TO 8°C.
IDEAL FOR SMALL SPACES.

FEATURES

Epoxy covered external case.
Interior ABS molded, resistant to chemicals with LED light built-in.
Internal fan for temperature homogenization.
Supports for adjustable in height shelves.

MODELO

Código	Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior)cm	Power W	Power HP	Nº Shelves	Weight Kg
2101286	105	146,5 28 27	182,5 36 38,5	110	1/8	8	65

Supplied with 8 steel trays, PVC laminated. For additional tray. Part No. **1001861**

CONTROL PANEL

Analog temperature controller.
On/Off switch internal light.
On/Off switch external light.

