



## Incubation chamber "PLX-IN"

FOR PLATELET BAGS "APL-54" SHAKING, CAPACITY 54 PLATELET BAGS.

FOR REGULABLE TEMPERATURES FROM 18 °C UP TO 60 °C.

RESOLUTION 0,1 °C. STABILITY  $\pm 0,1$  °C. HOMOGENEITY  $\pm 0,5$  °C. SET ERROR  $\pm 0,3$  °C.

***Peltier effect. No compressor group. Quiet, vibrations free, high accuracy and low consumption.***

### FEATURES

Made of transparent methacrylate with swinging door that allows an easy access to the working area and an inside chamber view during bags shaking.

AISI 304 stainless steel lower base.

Auto shut-off protective device in case of door opening and automatic restart in 30 seconds.

### CONTROL PANEL

Start-up main switch with digital temperature regulation.

### MODEL

Part No.	Height / Width / Depth (exterior) cm			Height / Width/ Depth (interior) cm			Power w	Weight Kg
<b>3000958</b>	44	63	65	33	60	50	360	38

It can be supplied for 110 V- 60 Hz on demand.



Note: For a correct bag preservation, it is recommended to keep a temperature between 20 and 24 °C.

*Incubation chamber "PLX-IN" Part No. 3000958  
with "APL-54" shaker Part No. 3000980*