



## Glass filtration equipment

### TECHNICAL SPECIFICATIONS

**Materials:**

100% borosilicate glass 3.3 except aluminium clamps.

**Capacities:**

Funnel 300 ml. Filtration flasks 1000 ml.

**Plate porosity:** Grade 2/3

**Sterilization:** Autoclave (max. 134°C), dry heat (max. 180°C).

**Diameter filters:** 47 mm

**Output connections to vacuum:** Tube 8 - 10 mm. diameter.

Complete equipment. **Part No. 5297023**

