

FEATURES

Compact and easy to use, designed for a maximum pressure load of 15 tons, especially for pressing in chemical research laboratories and materials for preparing ceramic composite. This machine has an integrated auto return hydraulic pump.

Tablets from 6.3mm to 19mm can be compacted with the MT matrix accessory, for KBR compression technique and fluorescence analysis works.

Cylinder elevation of 0-25mm.

MODEL

Part No.	Max. Pressure tons	Height / Width / Depth mm			Weight Kg	Max. Working area mm Height / Width / Depth
5810200	15	55	13	21	30	150 x 160 x 160

ACCESSORIES

- Part No. **5810201** MT Matrix 12.5mm
- Part No. **5810202** MT Matrix 6.3mm
- Part No. **5810203** MT Matrix 19mm
- Part No. **5810204** MT Matrix 3.2mm

