



Digital pH meter "pH-2004"

PORTABLE METER WITH IP65 PROTECTION
AGAINST WATER AND DUST.



FEATURES

Automatic temperature compensation and read out from 0 to 60 °C with temperature sensor Rt-2252. Results memory for up to 250 readings.

Protection against water and dust IP65 (Impermeable to water splashes from any direction.) Automatic buffer recognition and calibration for 1 or 2 buffers, recognition of up to 5 selectable buffers.

RS-232 output to connect a printer or a computer. Electrode connector: BNC type.

Includes: carry case, pH electrode and temperature sensor Rt-2252.

Direct connection to printer Part No. 4120113 (see page 247).

Universal arm support not included. Part No. 1001552 (see page 248).



Part No. 1001552

SPECIFICATIONS

	Range	Precision	Resolution
pH:	0 a 14	±0.02	0.01
mV:	±1800	±2	1
T ^a :	0 to 60	±0.5	0.1

Impedance: 1 x 10¹² Ohm.

MODEL

Part No.	Height / Width / Depth cm	Power source	Weight Kg
4120400	20 10 4.5	4 x 1.5V AA Batteries	0.5

ACCESSORY

pH probe and additional temperature sensor, gel filled. Temperature sensor Rt 2252.

With IP65 protective connectors.

Part No. **4120401**

