



# Digital pH meter "pH-2005"

**BENCH TOP MODEL WITH LARGE SCREEN.**

## FEATURES

Automatic or manual temperature compensation from 0 to 99.9 °C Automatic buffer recognition and calibration of 2 buffers, recognition of up to 3 buffers: pH 4, pH 7 and pH9.21.

RS-232 output to connect a printer or a computer.

Electrode connection BNC.

Temperature sensor connection Rt-2252.

Includes electrode's support stand.

## SPECIFICATIONS

	Range	Precision	Resolution
pH:	0 to 14	±0.02	0.01
mV:	±1999	±1	1
T <sup>a</sup> (°C):	0 to 99.9	±0.5	0.1
Impedance:	1 x 10 <sup>12</sup> Ohm.		



## MODEL

Part No.	Height / Width / Depth cm	Weight Kg
<b>4120500</b>	<b>30 20 7.2</b>	<b>1.5</b>

See accessories for electrode, temperature sensor and buffer solutions.  
Not compatible with printer Part No. 4120113 (see page 247).