

COMECTA Manual pipette dispenser "Pipet"

FOR ANY KIND OF PIPETTES FROM 1 TO 100 ML.

FEATURES

Ergonomic, comfortable, solid, lightweight and ease of use.

It has been designed with an ultracompressor silicone bulb that provides a soft manual control of the pipette, either in aspiration or dispensation.

Once the bulb is squeezed, control with the toggle both speeds, the aspiration and dispensation one. Push the middle button to expel blowing-out any liquid remaining in the pipette. 0,45 μm .

Measure: 25 mm \emptyset .

Weight: 130 g.

Manual pipette dispenser "Pipet". Part no **5900800**

COMPLEMENT

5 hydrophobe filters pack of 0,45 μm .

Part no: **5900801**

