



Automatic bag sealer “Seal-Com 2”

ACCORDING TO DIN 58953 SECTION 7 STANDARD. FOR POLYPROPYLENE-POLYESTER BAGS WITH PAPER BASE TO STERILIZE. CONTINUOUS SEALING. MANUAL OR AUTOMATIC ADVANCE.

Continuous Thermosealer

FEATURES

Exterior in enamel protected steel.
Heating time: 4'.
Sealing width: 10 mm.
Sealing speed: 6 m / min.
Acoustic signal of cycle starting.
Electronic regulation and digital reading of temperature.

CONTROL PANEL

Start-up light switch.
Manual or automatic start switch.

MODEL

Part No.	Height / Width / Depth (sealer) cm	Power W	Weight Kg
4100044	16 48 26	720	13

NEW
DESIGN

