



## Bag sealer “Seal-1000”

FOR POLYPROPYLENE AND POLYESTER BAGS WITH A STERILE PAPER BASE.

### FEATURES

Metal construction with stainless steel tray.

Pressure sealer unit with lateral lever with two bilateral cutting edges.

Constant sealing temperature (from 180 °C and 200 °C).

External support for rolls of up to 240mm wide, with adjustable separator sides to accommodate different roll sizes.

### CONTROL PANEL

Illuminated start switch.

Indicator, optimum seal temperature.

Part No.	Seal width mm	Maximum seal length mm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
<b>4101712</b>	<b>10</b>	<b>260</b>	<b>30 36 43</b>	<b>400</b>	<b>6.5</b>

