

Multifunctional grinding mill "TR-20"

FOR CRASHING AND SHEARING (INDUSTRIAL DESIGN)

APPLICATIONS

For a wide range of products grinding, solid, fragile, soft and semi-hard ones. For the use in pharmaceutical, chemistry and alimentary industry in general.

FEATURES

Made of stainless steel AISI 304.

Device for different kind of easily changeable sieves. From 0,12 mm Ø to 0,80 mm Ø.

It is supplied with a sieve of 0,18 mm Ø.

Only suitable for grinding samples with relative humidity below 20%.




Model	Part No.	Capacity Kg / h	Db	RPM	Voltage V	Measures (cm) Height / Width / Depth.	Power Kw	Weight Kg
TR-20	5810000	10 a 30	85	4200	380-50 Hz phase	70 40 80	2,2	85

ACCESSORIES

Sieves

Part no.	Sieve Ø mm.
5810007	1 mm
5810001	0,80 mm
5810002	0,38 mm
5810003	0,25 mm
5810004	0,18 mm
5810005	0,12 mm



SPARE PARTS Shear Part no. **5810006**

Device for an easy sieve insertion.