



Vacuum sealer "Sealcom-V"

OPTIMUM QUALITY SEAL AND PRESERVATION.
SUITABLE FOR SEALING SOLIDS AND LIQUIDS.

APPLICATIONS

Vacuum seal for any type of thermoplastic material: polyethylene, polypropylene, polyester, aluminium etc. Ideal for sealing solids, liquids, powders, granules, seeds, fruit, food, medical apparatus, forensics, chemical products, sludge samples, mud, cores, etc. Prevents the oxidation, mould and decomposition. Prolongs the quality, durability and freshness of items to be stored.

FEATURES

Bench top. Stainless steel AISI 304 robust construction with transparent lid. Easy to use and program. After closing the door the complete process is automatic. Operates using a oil vacuum pump (type AV-30 or N32).

CONTROL PANEL

Vacuum gauge with scale 0 to -0.1 Mpa.

Sealer temperature adjustment, high, medium or low.

Push button emergency stop.

LED indicator of operation.

Digital selection and display of the following parameters:

- Vacuum time 0-99 seconds.
- Seal time 0 to 9.9 seconds
- Discharge time 0 to 9.9 seconds



| MODEL | Part No. | Seal Width mm | Max. seal length mm | Pump rate m ³ /h | Vacuum limit Mpa | Height / Width / Depth (usable) cm | Height / Width / Depth (exterior) cm | Power W | Weight Kg |
|-----------|----------|---------------|---------------------|-----------------------------|------------------|------------------------------------|--------------------------------------|---------|-----------|
| SEALCOM-V | 4100050 | 10 | 260 | 10 | -0.1 | 5 28 39 | 36 35 52 | 370 | 35 |

Supplied with a maintenance kit that includes 2 heaters, 3 PTFE heater covers. 1 lid seal and a 250 ml bottle of oil.

ADDITIONAL SPARES

Maintenance kit for Sealcom-V: 2 heaters, 3 PTFE heater covers. 1 lid seal. Part No. 4100051
Vacuum pump oil: 1 litre bottle (Type AV-30). Part No. 4100052