



## Cold extraction unit "EF-1425"

### DESCRIPTION

Specifically designed unit for cold extraction of up to 6 samples simultaneously, using acetone or other solvents. Complements the Dosi-Fiber for the analysis of crude fibre.

Components: 1 litre receptacle for solvent recovery, suction tube and 6 adapters to hold the filtration crucibles Part No. 4000601.

A pump or jet pump will be required to connect to the suction tube.

Comes complete with: blank stoppers in the case that not all of the 6 crucible places are used and a 1 litre solvent recovery vessel.

### FEATURES

Stainless steel AISI 304, extractor, modular construction, easy to use and fast removal for efficient cleaning. Support and crucible rack, made from hardened glass.

### MODEL

Part No.	Height / Width / Depth (exterior) cm	Weight Kg
<b>4001425</b>	<b>45 20 34</b>	<b>3</b>

### ACCESSORY

Jet pump.

Part No. **7000293**

