

## APPLICATIONS

Microscopy illumination for small samples free from colour alteration.

Especially recommended for microphotography of live biological samples, cultures, micro-organisms, process control, electronic components, dental studies, etc.

## FEATURES

Halogen lamp 150 W. Adjustable illumination intensity control. Forced air circulation. Heat absorption filter.

Safety transformer fitted. Primary fuse protection.

## CONTROL PANEL

Mains switch. Potentiometer controlled light intensity.

## MODEL

Part No.	Height / Width / Depth exterior cm	Power W	Weight Kg
<b>5901016</b>	<b>15 19 25</b>	<b>150</b>	<b>4</b>

**SPARE Lamp 150 W.**

Part No. **5901017**

Supplied with double flexible arm 50 cm long.

