



Measuring range “COM-PIP”

FEATURES

Made of polypropylene with silicon joints, for volumes up to 25 ml. Precise and safe, ergonomic design that allows filling and emptying without having to change the setting. Flexible thermoplastic seal against any type of burette. A button on the side empties the burette automatically. Easily dismantled and cleaned.

Part No.

1001131

Capacity

0 to 2 ml

Identification colour

Blue

1001132

0 to 10 ml

Green

1001133

0 to 25 ml

Red



PIPETTE RACK “COM-PIP”

Made from AISI 304 stainless steel.

Capacity 6 pipettes of differing sizes.

Part No. **1001702**

