



Colony counter "Digital S"

FOR PETRI DISH OF UP TO 100 mm DIAMETER.
LED ILLUMINATION.

APPLICATIONS

Molecular Biology, Genetics and microbiology.

FEATURES

The unit records each pulsed colony count on the digital screen to a total of 999 counts.

Pressure to count can be made by any way (pencil, or pen tip etc.).

Epoxy coated casing. Reversible interior background part in white or black.

Steady, (glare free) LED illumination, clear and wide view visibility.

Removable plates for housing varying sizes of Petri dishes.

Location for reversible contrast plates.

Magnifying glass X2 with adjustable height and rotation support

CONTROL PANEL

Mains switch.

Numerical digital display: 3 digits.

Audible confirmation beep for each count.

Zero reset button. RESET.

MODEL

Part No.	Ø of magnifying glass cm	Magnification	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4905002	12.5	2x	25 23 30	20	3.5

Supplied complete with Wolffhuegel disc.



ACCESSORY

5 Wolffhuegel discs.

Part No. **4905001**