



Heating mantle for flasks

ADAPTABLE FOR THE "AGIMAN" STIRRER UNIT.

MODELS

Part No.	Flask capacity ml	Ø / Height (Exterior dims.) cm		Power W	Weight Kg
3000447	250	15	6.5	130	0.8
3000448	500	17	7	270	0.9
3000449	1000	20	8.5	410	1
3000450	2000	24	11	530	1.5

ACCESSORIES

Mantle magnetic stirrer "Agiman". With interchangeable mantle adapters. Electronic control of stir speed (up to 1600 r.p.m.). Part No. **7000446** See page 26. With connection at the back for an electronic contact thermometer "Sensoterm" or digital electronic regulator "Electemp-TFT". See controller (pages 293 and 294).



Stirring mantle with adapter device.



Stirring mantle, complete.