



# Stackable Anatomy - Histology sample storage chest

## FEATURES

Manufactured in steel epoxy coated. Resistant to corrosive atmospheres.

## STACKABLE MODULAR SYSTEM

### Cabinet for storing slides.

Complete with 14 drawers for slides of 26 x 76 mm.

Drawer capacity: 350 slides.

Part No. **1052014**

### Cabinet for storing transparencies.

Complete with 7 drawers for transparencies of 50 x 60 mm.

Drawer capacity: 230 transparencies.

Part No. **1005207**

### Cabinet for storing paraffin wax mounted samples.

Complete with 10 drawers.

Part No. **1052010**

## MODELS

Part No.	Height / Width / Depth (usable) cm			Height / Width / Depth (exterior) cm			Weight Kg
<b>1052014</b>	9	2.6	45	14	48	48	9.5
<b>1005207</b>	7.2	5.8	47	12	48	48	6.5
<b>1052010</b>	2.6	43.5	47	34.5	48	48	18.5

### Base and lid units

Mounted on the top or as a support base. Fits the individual modules to form distinct user definable cabinets.

Part No.		Height / Width / Depth (external) cm			Weight Kg
<b>1005208</b>	<b>Base</b>	5,5	48	48	2,5
<b>1005209</b>	<b>Lid</b>	5,5	48	48	2,5

