

# Multichannel variable volume dispensing pipettes “Multi-Var 8”

NUMBER OF CHANNELS: 8.

## FEATURES

Due to its excellent tips alignment, they offer a great advantage in microtitrer techniques such as: agglutination, immunological radio assays, hemagglutination, antibiotic test, etc. (All models are calibrated according to DIN 12600 recommendations)

MODELS	Part No.	Tips Range µl	Nominal Volume Accuracy		Coefficient variation Nominal volume		Increments µl	Tips types µl
			A < ±%		CV < %			
			Min %	Max %	Min %	Max %		
Multi-Var 8	5905027	5-50	0.8	6	0.4	3	0.1	200
Multi-Var 8	5905028	10-100	0.8	4	0.3	2	0.2	200/300
Multi-Var 8	5905029	20-200	0.8	4	0.3	1.5	0.2	200/300

## ACCESSORIES

Pipette tips rack. Part No. **5490037**

Carousel support for a maximum 3 multichannel pipettes. Part No. **5490038**

